



Micrometrix™



Model PCA PARTICLE CHARGE ANALYZER

Features

- Digital Display
- LED Indicators
- milli-Volt Output
- Zero Offset
- Sensitivity Adjustment
- Portable

Charge Measurement and Charge Demand Analysis



The Micrometrix™ Particle Charge Analyzer measures ionic and colloidal charge in liquid samples. Optimum chemical dosage can be determined much quicker than with standard “jar testing”. Anionic or cationic Charge demand is determined by titration. The digital display indicates the zero charge endpoint value when titrating.

Applications

- Water Treatment
- Paper Industry
- Wastewater
- Chemical Industry
- Ceramics
- Minerals
- Remediation

Benefits

- Establish Chemical Dose
- Optimize Process Performance
- Quality Control
- Replace Jar Tests
- Minimize Chemical Cost
- Reduce Residuals
- Quick and Easy to Use



770.271.1330

micrometrix.com

contact@micrometrix.com


©Micrometrix Corp 2012

Engineering Specifications

The Streaming Current cell (SC cell) determines the charge of the sample and the end point of the titration. The measurement cell consists of a precision bore cylinder closed at the bottom end and containing two electrodes, one at the bottom, and an upper electrode. The electrodes are connected to the contacts extending from the lower front portion of the instrument housing. The measurement cell is designed as a container to allow sample to be poured in from the top. The typical sample volume is 100ml. A precision piston oscillates up and down in the cylinder with a frequency of approximately 4hz. Polyelectrolyte (polymer) or coagulants are used as titrants to determine the charge demand of the sample. Colloidal particles are temporarily attached to the piston and cylinder walls. The mobile counterions of the fixed electrolyte move through the liquid stream creating an electric current due to the partial charge distribution measured between the two electrodes. This streaming current is measured by the electronics in the main unit.

The instrument shall be a Particle Charge Analyzer Model PCA manufactured by Micrometrix Corporation.

| Specifications | Model PCA |
|------------------------------|---|
| Measurement | Streaming Current |
| Power | 110 Vac / 230 Vac (Optional) |
| Range | -1000 to +1000 mV |
| Accuracy | 0.1% |
| Display Type | LED - Digital |
| Sample Size | 10 ~ 100 ml |
| Connection Type | Piston connected magnetically |
| Response Time | 1 Second |
| Self Diagnostic | LED |
| Electrode and Cell Materials | 316 SS, PTFE |
| Outputs | -1000 to +1000 mV |
| Options | Endpoint Titrator, Carrying Case, Data Logger |
| LED Indicators | Anionic and Cationic |
| Zero Adjust | Full Range |
| Type / Use | Portable/ Benchtop |
| Dimensions | 7" x 7" x 14" |
| Weight | 20 lbs , 9 kg |



770.271.1330

micrometrix.com

contact@micrometrix.com

©Micrometrix Corp 2012